## **CONTEC**











CMS60D-VET Pulse Oximeter



#### SpO<sub>2</sub>

Measurement Range: 0~100%

Resolution: 1%

Accuracy: 70~100%, ±2%; Below 70%, unspecified

Alarm Range: 0~100%

#### Pulse Rate

Measurement Range: 30~250bpm

Resolution: 1bpm

Accuracy: ±2 bpm or±2% (select larger)

Alarm Range: 0~254bpm

#### Storage Environment

a) Temperature: -40°C ~ +60°C

b) Relative humidity: ≤95%

c) Atmospheric pressure: 500hPa ~ 1060hPa

### Operating Environment

a) Temperature: 10°C ~ +40°Cb) Relative Humidity: ≤75%

c) Atmospheric pressure: 700hPa ~ 1060hPa

#### Feature

Display Mode: 1.8" OLED display

Power Supply: Two 1.5V (AA size)alkaline battery

Battery Voltage: DC 2.6V~3.6V

SpO2, Pulse Rate, Bar Graph, Pulse Waveform display

Battery capacity indication

Brightness of screen can be changed

Pulse sound indication

# CMS60D-VET Pulse Oximeter

- With alarm function
- Synchronous Display on the Oximeter and Computer used USB
- Inner flash memery for 24 hours data storage, and the storage data can be uploaded to computers
- With review function
- Clock function
- It can be connected with an external oximeter probe
- Small in volume, light in weight and convenient in carrying
- Application: The device is suitable for being used in family, pet hospital. (It is recommended to use the device when the animal is still.)

**Илан Медицинска Апаратура ООД** гр. Варна, ул. Кирил Шиваров 9 Б тел. **0700 17373** факс 052 612258 e-mail: office@ilan.bg