



reddot design award  
winner 2011

## Real Time Interaction



### Instant Examination View



# On-the-spot examination

## Power / Freeze image on one single button

- Simply press the indicated button for 2 seconds to turn ON and OFF the device
- Quick press the indicated button to capture and release the display image, when the device is ON.



### **MDscope meets with safety standards include:**

Medical Electronic Equipment EN60601-1-2  
Medical Device Directive 93/42/EEC  
FDA clearance  
CE certificate

### **Product Specification**

Body size: 2.5" width x 8" height  
2.4" color LCD display  
Light weight: 218 grams (0.48 lb)  
Battery lifespan: 4 hours (AA x 2)

### **Special features include:**

Image Freeze function  
Video out port

APPLE BioMedical Inc.  
MDSCOPE REF: 101  
Made in Taiwan  
ELECTRONIC  
FSM 12071458  
CE  
FCC  
RoHS



**Detachable camera probes for different examinations as needed.**



Patents pending worldwide


**Video port to provide simultaneous image to TV/ monitor, or to download image to Mac / PC system.**



*monitor / TV*



*MDscope image dongle  
for MAC / PC system*



## Packages



### MS101 Elite Pack

Include 1 standard 45mm camera probe



### MS101 Deluxe Pack D1

Include 3 camera probes – 45mm / 75mm / 150mm  
And 1 image dongle with USB cable

### MS101 Deluxe Pack D2

Include 3 camera probes – 45mm / 75mm / 300mm  
And 1 image dongle with USB cable

## Accessories



### MS101-002

Standard camera probe 45mm  
For ear examination.



### MS101-012

Long camera probe 75mm  
For nose examination.



### MS101-013

Flexible camera probe 150mm  
For throat examination.



### MS101-015

Disposable specula 5mm tip  
250 pcs per box.



### MS101-011

Image dongle with USB cable,  
Compatible with Mac and PC.



### MS101-014

Flexible camera probe 300mm  
For long distance examination.



# ИЛАН Медицинска Апаратура ООД

централен офис: Варна 9000, ул. Кирил Шиваров 9 Б  
тел. 0700 17373, тел./факс 052 612258, GSM 0898 587197, 0885 565080  
<http://www.ilan.bg>, e-mail: [office@ilan.bg](mailto:office@ilan.bg)